|                |          |   |             |
|----------------|----------|---|-------------|
| <b>STERILE</b> | <b>A</b> | Sterilized using aseptic processing techniques (filtration) | <b>CE</b>   |
|                |          |   | <b>2797</b> |

**Material included:**

***GM501 Gradient 45 %***

Product code: 4 GM 501G-45-10

- 1 x 10 ml GM501 Gradient 45 %

Product code: 4 GM 501G-45-50

- 1 x 50 ml GM501 Gradient 45 %

Product code: 4 GM 501G-45-100

- 1 x 100 ml GM501 Gradient 45 %

Product code: 4 GM 501G-45-250

- 1 x 250 ml GM501 Gradient 45 %

***GM501 Gradient 90 %***

Product code: 4 GM 501G-90-10

- 1 x 10 ml GM501 Gradient 90 %

Product code: 4 GM 501G-90-50

- 1 x 50 ml GM501 Gradient 90 %

Product code: 4 GM 501G-90-100

- 1 x 100 ml GM501 Gradient 90 %

Product code: 4 GM 501G-90-250

- 1 x 250 ml GM501 Gradient 90 %

**Material needed but not provided:**

- 3cc syringes with 1 1/2" 21 g needle

**Product specifications and quality control:**

- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for GM501 Gradient media are available upon request and can also be downloaded from our website (www.gynemed.de).
- Not MEA tested
- GM501 Gradients are manufactured and tested according to the following specifications:***
  - Chemical composition
  - pH release criteria: 7.20-7.60, pH stability criteria: 7.20-7.90
  - Sterility test by the current Ph. Eur. 2.6.1 / USP <71>; *No growth*
  - Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): *< 0.5 EU/ml*
  - Density: *1.0500-1.0700 g/ml (Gradient 45%)*  
*1.1050-1.1150 g/ml (Gradient 90%)*
  - Osmolality: *310-340 mOsm/kg (Gradient 45%)*  
*320-350 mOsm/kg (Gradient 90%)*
  - Human Sperm Survival Assay: *≥ 80% survival after 4 hours exposure of density selected spermatozoa to the test medium; ≥ 75% survival after 24 hours exposure of density selected spermatozoa to the test medium*

**Calculation of G-forces:**

- The G-force of your centrifuge can be calculated using this formula:

g = 1.118 x r x rpm² or rpm = Square root (g/ (1.118 x r))

r = radius of centrifuge in mm

rpm = rotations per minute/1000

**Example 1**

r = 150 mm

rpm = 1200 rotations per minute

g = 1.118 x 150 x 1.44 = 242 g

**Example 2**

r = 150 mm

g = 300 g

rpm = SQR (300/ (1.118 x 150)) = 1.33

rpm = 1330 rotations per minute

**Used abbreviations:**

- EBSS** Earle’s Balanced Salt Solution
- ICSI** IntraCytoplasmatic Sperm Injection
- IVF** In Vitro Fertilization
- IUI** Intra Uterine Insemination

|  |  |  |  |  |  |
|--|--|--|--|--|--|
|  |  |  |  |  |  |
|--|--|--|--|--|--|

**EN** Intended use/Intended users:

- GM501 Gradient 45/90 is a ready-to-use gradient system for semen prepara-tion. The following gradients are available:
  - Gradient 45%
  - Gradient 90%
- GM501 Gradients can be used in combination with IUI, IVF and ICSI.
- For professional use only.

**Composition:**

- GM501 Gradients consist of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle’s balanced salt solution).

**Warning before use:**

- Do not use the product if it becomes discolored (medium contains phenol red), cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.
- Do not use the product if the expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.

**Instructions for use:**

- Importantly: each laboratory should consult its own optimized and validated laboratory procedures.

**Instruction for use with fresh semen samples:**

Mix the density gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37°C
- Transfer 2.5 ml of Gradient 45% into a sterile disposable centrifuge tube.
- Using a 3cc syringe with a 1 1/2" 21 g needle, place 2.5 ml of Gradient 90% under the Gradient 45%. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle on the bottom of the test tube and slowly dispensing the Gradient 90%. This two-layer gradient is stable for up to two hours.
- Gently place up to 2.5 ml of liquefied semen onto the Gradient 45% using a transfer pipette or syringe.
- Centrifuge for 15 to 18 minutes at 350 g to 400 g. When this centrifugation is completed, you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
- Remove supernatant down to the pellet.
- Using a syringe, add 2-3 ml of sperm washing medium and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300 g.
- Higher sperm concentration will require the maximum 10 minutes centrifuga-tion to ensure a complete and thorough sperm wash.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate me-dium.

**Instructions for use with frozen semen samples:**

Mix the density gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37°C.

- Transfer 1 ml of Gradient 45% into a sterile disposable centrifuge tube.
- Using a 3cc syringe with a 1 1/2" 21 g needle, place 1 ml of Gradient 90% under the Gradient 45%. Take care that the two layers are distinctly separated. This is done by placing the tip of the needle on the bottom of the test tube and slowly dispensing the Gradient 90%. This two-layer gradient is stable for up to two hours.
- Gently place the thawed semen sample onto the Gradient 45% using a trans-fer pipette or syringe (0.5 ml maximum).
- Centrifuge for 15-20 minutes at 350 g.
- Remove supernatant down to (but not below) the 0.5 ml mark above the pellet.
- Using a syringe, add 2-3 ml of sperm washing medium and resuspend the pellet.
- Centrifuge for 8 to 10 minutes at 300 g.
- Remove supernatant down to the pellet and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate me-dium.
- Centrifuge for 8 to 10 minutes at 300 g.
- Remove supernatant and repeat steps 7 and 8.
- Remove supernatant and replace with a suitable volume of appropriate me-dium.
- Keep away from (sun)light

*If samples do not liquefy and therefore do not pass through the layers, increasing the centrifugal force up to, but no more than, 500 g will help to separate the sperm.*

**Storage and disposal instructions:**

- The shelf life is 18 months from time of manufacture.
- Store products between 2-25°C before first use, once opened at 2-8°C.
- The product can be used safely up to 7 days after opening, if sterile conditions are maintained, and the products are stored at 2-8°C.
- Stable after transport at elevated temperatures (up to 5 days at ≤ 37°C)
- The devices need to be disposed in accordance with local regulations for dispo-sal of medical devices.
- Keep away from (sun)light

**Precautions and warnings:**

- All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic technique must be used to avoid possible contamination.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred must be reported to Gynemed and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.
- Only for the intended use.

**Clinical benefit:**

- The clinical benefit of Density Gradient Media is related to the fact that these media are used during ART-procedures of patients and couples with infertility problems.
- It is expected that the outcomes of ART procedures in which density gradient is used are comparable with the ART outcomes (i.e. pregnancy rate and deli-very rate) as published in the annual peer-reviewed report of the European Society of Human Reproduction and Embryology (ESHRE).
- For further questions regarding the safety and performance, please contact Gynemed GmbH & Co. KG for customer or technical support.

**Safety and clinical performance**

- For questions regarding the safety and performance, please contact our Gy-nemed regulatory affairs or technical support team.

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| Symbols glossary                  |   |  |   |
|-----------------------------------|---|--|---|
| Symbols as defined in ISO 15223-1 |   |  |   |
| <span></span>                     |   | <span></span>  |   |
| <b>REF</b>                        | Catalogue number  | <b>LOT</b>   | Batch code  |
| <span></span>                     |   | <span></span>  |   |
| <span></span>                     | Use-by date   | <span></span>  | To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. |
| <span></span>                     | Temperature limit   | <b>STERILE A</b>   | Sterilized using aseptic processing techniques  |
| <span></span>                     |   | <span></span>  |   |
| <span></span>                     | Keep away from (sun)light   | <span></span>  | Do not resterilize  |
| <span></span>                     |   | <span></span>  |   |
| <span></span>                     | Consult instructions for use  | <span></span>  | Do not use if package is damaged or opened and the user should consult the Instructions for Use for additional information.   |
| <span></span>                     |   | <span></span>  |   |
| <span></span>                     | Date of manufacture   | <span></span>  | Sterile barrier system.   |
| <span></span>                     |   | <span></span>  |   |
| <span></span>                     | Manufacturer  | <span></span>  | Single sterile barrier system with outer protective packaging   |
| <span></span>                     |   | <span></span>  |   |
| <b>MD</b>                         | Medical Device  | <b>UDI</b>   | Unique device identifier  |
| <span></span>                     |   | <span></span>  |   |
| Symbol as defined in MDR 2017/745 |   |  |   |
| <span></span>                     |   | <span></span>  |   |
| <b>CE</b><br>2797                 | CE marking by Notified Body 2797  | <span></span>  |   |
| For technical support             |   |  |   |
| <span></span>                     | GYNEMED GmbH & Co. KG<br>Wagrienring 24b<br>23730 Sierksdorf<br>- Germany - | Phone: + 49 (0) 4561 5 13 25-00<br>Fax: + 49 (0) 4561 5 13 25-19 | E-mail: <a href="mailto:info@gynemed.de">info@gynemed.de</a><br>URL: <a href="http://www.gynemed.de">www.gynemed.de</a>   |

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