**Material included in the package:**

- Product code: 4 HSA-0005
- 1 x 5 ml GM501 HSA (10% human serum albumin)

**Material not included in the package:**

- CO<sub>2</sub>-incubator at 37°C
- LAF bench (ISO 5 environment)

**Product specifications and quality control:**

- GM501 HSA is obtained from blood of donors of certified blood banks according to international guidelines. The HSA is dialyzed and subsequently solved in saline buffer.
- The protein supplement is manufactured using the modified Coon-Ohey cocktail ethanol fractionation procedure and embryo manipulation with filtration and additionally heated for 10-11 hours at 60°C.
- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respect to lot number.
- The MSDS for GM501 HSA is available upon request and can also be downloaded from our website.
- GM501 HSA is manufactured and tested according to the following specifications:

pH (at 37°C, 6% CO <sub>2</sub> )	7.20-7.40
Osmolality (mOsm/kg)	265-285
Sterility	sterile - SAL10 <sup>-3</sup> (Sterility Assurance Level)
Endotoxins (EU/ml)	< 1.00
MEA	≥80
(1-cell assay, blastocysts after 96 h in %)	

.....

**EN: Intended use/intended users:**

- GM501 HSA (10 % human serum albumin) is intended for use in assisted reproductive procedures involving gamete and embryo manipulation. These procedures include the use of GM501 HSA solution as a supplement for culture medium.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors)

**Compositio:**

- 10 % (100mg/ml) GM501 HSA in saline buffer.

**Instructions for use:**

- To supplement the 50 ml culture medium (e.g. GM501 Basic or GM501 Basic with Gentamicin) with GM501 HSA to a protein end concentration of 1 % .Keep 45 ml culture medium in the bottle and add 5 ml of GM501 HSA. Mix well.

**Precautions and warnings:**

- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured to European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAF-bench) to avoid possible contamination.
- Only use the intended use as described in the instructions for use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

**Pre-use checks:**

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective when the product is delivered.

**Storage instructions and stability:**

- The shelf life is 12 months from time of manufacture.
- Store between 2-8 °C.
- Protect from light.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Stable after transport (maks. 5 days) at elevated temperature (≤ 37°C).
- Do not use after expiry date.

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GYNEMED GmbH & Co. KG