GM50 Wash

Material included in the package:

- 4 GM 501W-20
- 1 x 20 ml GM501 Wash
- Product code: 4 GM 501W-50
- 1 x 50 ml GM501 Wash
- Product code: 4 GM 501W-500
- 1 x 500 ml GM501 Wash

Material included in the package:

- GM501 Wash with Phenolred and Gentamicin
- Product code: 4 GM 501W+G+PR-20
- 1 x 20 ml GM501 Wash with Phenolred and Gentamicin
- Product code: 4 GM 501W+G+PR-50
- 1 x 50 ml GM501 Wash with Phenolred and Gentamicin

Material not included in the package:

- CO₂-incubator at 37°C
- Plain dishes
- LAf bench (ISO 5 environment)
- Microscope
- Test tubes

Product specifications and quality control:

- All raw materials are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
- A certificate of analysis is available for each batch upon request from our website with respective lot number.
- The MSDS for the media are available upon request and can also be downloaded from our website.
- The media are manufactured and tested according to the following specifications:

pH (at 37°C, 6% CO ₂)	7.20-7.50
Osmolality (mOsm/kg)	270-290
Sterility	sterile - SAL101 (Sterility Assurance Level)
Endotoxins (EU/ml)	<=0.25
Mouse Embryo Assay (MEA)	>20
(1-cell assay, blastocysts after 96 h in %, 60min exposure at zygote stage)	

EN: Intended use/Intended users:

- The media are ready-to-use media designed for washing procedures of human oocytes and embryos and any short-term handling procedures outside the incubator like washing after Hyaluronidase treatment (denudation), ICSI, polar body or blastomere biopsy and others.
- The intended users are IVF professionals (lab technicians, embryologists or medical doctors).

Composition:

GM501 Wash

NaCl, KCl, KH₂PO₄, MgSO₄·7H₂O, NaHCO₃, CaCl₂·2H₂O, D(+)-Glucose anhydrous, Na-Lactate (50 % solution), Na-Pyruvate, EDTA, Alaninyl-Glutamine, water, non-essential and essential Amino Acids, 15 mM HEPES, human serum albumin (5.00 g/liter).

GM501 Wash with Phenolred and Gentamicin

NaCl, KCl, KH₂PO₄, MgSO₄·7H₂O, NaHCO₃, CaCl₂·2H₂O, D(+)-Glucose anhydrous, Na-Lactate (50 % solution), Na-Pyruvate, EDTA, Alaninyl-Glutamine, water, non-essential and essential Amino Acids, 15 mM HEPES, human serum albumin (5.00 g/liter), Gentamicin (10 mg/liter), Phenolred.

Instructions for use:

- The medium must be equilibrated over night in a humidified CO₂-incubator (at 5-7 % CO₂, 37°C).

Precautions and warnings:

- GM501 Wash with Phenolred and Gentamicin contains the antibiotic: Gentamicin Sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot totally be excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin, manufactured in European Pharmacopoeia specifications by established processes.
- Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
- Always work under strict hygienic conditions (ISO 5 environment, e.g. LAf-bench) to avoid bacterial contamination, even when GM501 Wash medium contains Gentamicin.
- Only for the intended use.
- The user facility of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability, where applicable.

Pre-use checks:

- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
- Do not use the product if seal of the container is opened or defective

Storage instructions and stability:

- The shelf life is 6 months from time of manufacture.
- Store between 2-8°C.
- Do not freeze above zero.
- Keep away from (sun) light.
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at the elevated temperature (≤ 37°C).
- Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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