
# Vitrification Cryotop Method for Oocyte and Embryo Vitrification Media

Code:VT801

Symbols on the Label

June, 1, 2015 Ver.5

	Sterilized Using Aseptic Processing Techniques
	Sterilized Using Irradiation
	Do not reuse
	Catalog Number
	Batch Code
	Use by
	Storage Temperature
	CE mark
	Manufacturer
	Do not re-sterilize
	Do not use if package is damaged
	Consult Instructions for use

**EN: Intended use**

- Vitrification Media is to be used for vitrifying Oocytes and Embryos.
- CAUTION:** Sterile if the package is unopened or undamaged. Do not use if packaging is broken.

**Vitrification Media**

- No.0 Basic Solution (BS): 1 X 1.5ml vial (only for Oocyte Vitrification)
- No.1 Equilibration Solution (ES): 1 X 1.5ml vial
- No.2 Vitrification Solution (VS): 2 X 1.5ml vials

**Recommended equipment**

- Cryotop: 1 Cryotop stores up to 4 Oocytes or 4 Embryos as a recommendation.
- Repro Plate: with 6 wells

**Instructions for use (IFU)**

- Preparation**
- Fill 90% of the Cooling Rack with fresh liquid nitrogen.
  - Compare the width of perivitelline space with the thickness of zona pellucida and record it.
- NOTE:** Use a pasteur pipette that has a suitable internal diameter for Oocyte (External diameter: 120µm) or Embryo.

**EQUILIBRATION**

- CAUTION:** Equilibration procedures for Oocyte and Embryo are different.

**Oocyte Equilibration**

- Drop 20µl for BS and 300µl each for VS1 and VS2 on the Repro Plate with micro pipette.
- Transfer the Oocyte from the culture dish to the BOTTOM of BS.
- Step 1. Add ES 20µl gently to the TOP of BS and leave it for 3 minutes.
- Step 2. Add another ES 20µl gently to the TOP of BS and leave it for 3 minutes.
- Step 3. Add another ES 240µl gently to the TOP of BS and leave it for 9 minutes.

**Embryo Equilibration**

- Drop 300µl for ES, VS1 and VS2 on the Repro plate using micro pipette.
- Transfer the Embryo to the TOP centre of ES.
- It will spontaneously begin to shrink and then gradually return to its original size by absorbing the ES solution (within 15 minutes).

**VITRIFICATION**

- CAUTION:** The following steps from 4 to 9 should be completed between 60 and 90 seconds.
- Aspirate the Oocyte (Embryo) in ES at the tip of pasteur pipette. Transfer the Oocyte (Embryo) to the TOP center of VS1.
- Aspirate the Oocyte (Embryo) with pasteur pipette and blow it out. Repeat this process 3 times, changing the position in VS1.
- Transfer the Oocyte (Embryo) to VS2. Move the Oocyte (Embryo) changing positions twice in VS2 with pasteur pipette.
- Place the Oocyte (Embryo) by the black line on the Cryotop. Make a planar droplet.
- Check if the Oocyte (Embryo) is on the Cryotop with a minimal volume of VS2 (less than 0.1µl) under the microscope.
- Quickly plunge the Cryotop into the liquid nitrogen.
- Put the Cryotop in a cane and store it in a storage tank.

**Quality Control Specification**

Each lot of Vitrification Media receives the following tests:

- Sterility by the current USP Sterility Test <71> (Solutions)
- Endotoxin by LAL methodology
- Mouse Embryo Assay (One Cell)
- pH (Ph.Eur., USP)
- Osmolality (Ph.Eur., USP)

**Storage instructions and stability**

- Solutions: Store the vials at 2-8 °C
- The product is stable until the expiry date shown on the vial label.

**Composition**

- HEPES within Basic Culture Medium
- Ethylene Glycol
- Dimethyl Sulfoxide
- Trehalose
- Hydroxypropyl Cellulose

**Warning**

- Do not re-sterilize
- Do not use solution that shows cloudiness or has turned yellow.
- Do not use if sterile packaging is broken.
- Upon delivery media must be stored in original unopened container and refrigerated at 2-8 °C.
- To avoid contamination, do not reuse.

**Cautions**

- Read the instructions for use prior to use.
- This product is intended to be used by medical specialists trained in fertility treatment.
- Aseptic technique should be used.
- Use sterilized equipment and materials only.
- Decontaminate the workstation.
- Follow procedures in an environmentally controlled room.
- The long term safety is unknown.

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**References**


- Cobo A., Obstetric and perinatal outcome of babies born from vitrified oocytes. Fertility and Sterility, 2014.
- Rienzi L., Consistent and predictable delivery rates after oocyte vitrification: an observational longitudinal cohort multicentric study. Human Reproduction, 2012.
- Shi W., Perinatal and neonatal outcomes of 494 babies delivered from 972 vitrified embryo transfers. Fertility and Sterility, 2012.
- Cobo A., Outcomes of vitrified early cleavage-stage and blastocyst-stage embryos in a cryopreservation program: evaluation of 3,150 warming cycles. Fertility and Sterility, 2012.
- Trokudes KM., Comparison outcome of fresh and vitrified donor oocytes in an egg-sharing donation program. Fertility & Sterility, 2011.
- Inoue F., Hydroxypropyl cellulose as a macromolecular supplement for cryopreservation by vitrification of bovine oocytes and blastocysts and human oocytes. ESHRE and ASRM, 2011.

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CE 0086

## Code:VT801

STERILE A	Sterilized Using Aseptic Processing Techniques
STERILE R	Sterilized Using Irradiation
Do not reuse	
REF	Catalog Number
LOT	Batch Code
Use by	
Storage Temperature	
CE 0086	CE mark
Manufacturer	
Do not resterilize	
Do not use if package is damaged	
Consult Instructions for use	

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Copyright©KITAZATO CORPORATION All Rights Reserved.
Copyright © KITAZATO CORPORATION All Rights Reserved.
KITAZATO BioPharma Co., Ltd.