

Vitrification Cryotop Method for Oocyte and Embryo Thawing Media

CE 0086

Code:VT802

TS Warming
(+42000°C/min)
1min, 37°C

Symbols on the Label

June, 1, 2015, Ver.5

	Sterilized Using Aseptic Processing Techniques
	Sterilized Using Irradiation
	Do not reuse
	Catalog Number
	Batch Code
	Use by
	Storage Temperature
	CE mark
	Manufacturer
	Do not resterilize
	Do not use if package is damaged
	Consult Instructions for use

EN: Intended Use

- Thawing Media are for thawing vitrified Oocytes and Embryos.
- CAUTION:** Sterile if the package is unopened or undamaged. Do not use if package is broken.

Thawing Media

No.1 Thawing Solution (TS): 2 x 4.0ml vial
No.2 Diluent Solution (DS): 1 x 4.0ml vial
No.3 Washing Solution (WS): 1 x 4.0ml vial

Recommended equipment
Petri Dish: 35mm for TS
Repro Plate with 6 wells

Instructions for Use

Preparation

- Warm TS vial with a cap and a Petri Dish in an incubator to 37°C. Pour the full content of TS into the Petri Dish.
- Drop 300µl each of DS, WS1 and WS2 on the Repro Plate with micro pipette.

NOTE: Use Pasteur pipette that has a suitable internal diameter for Oocyte (External diameter: 120µm) or Embryo.

Thawing

- Quickly immerse the Cryotop strip completely into the TS. Leave it for 1 minute.
- Aspirate the Oocyte (Embryo) with the Pasteur pipette and gently place it on the BOTTOM of the DS. Leave it for 3 minutes.
- Aspirate the Oocyte (Embryo) with the Pasteur pipette and gently place it on the BOTTOM of the WS1. Leave it for 5 minutes.
- Aspirate the Oocyte (Embryo) with the Pasteur pipette and gently place it on the TOP of the WS2. After the Oocyte (Embryo) drops to the bottom of the WS2, repeat this process twice.
- Transfer the Oocyte (Embryo) to a culture dish containing the appropriate culture medium. Incubate the Oocyte (Embryo) in a 37°C incubator to complete recovery.

NOTE: Oocytes for 2 hours and Embryos for 3 hours as an idea.

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Quality Control Specification

Each LOT of Thawing Media receives the following tests:

- Sterility by the current USP Sterility Test <71> (Solutions)
- Endotoxin by LAL methodology
- Mouse Embryo Assay (One Cell)
- pH (Ph. Eur., USP)
- Osmolarity (Ph. Eur., USP)

Storage Storage instructions and stability

Solutions: Store the vials at 2-8°C

The products are stable until the expire dates shown on the labels.

Composition

- HEPES within Basic Culture Medium
- Trehalose
- Hydroxypropyl Cellulose

Warning

- Do not re-sterilize
- Do not use solution that shows signs of cloudiness or has turned yellow.
- Do not use if sterile packaging is broken.
- Upon delivery media must be stored in original unopened container, refrigerated at 2-8°C.
- To avoid contamination, do not reuse.

Cautions

- Read the instructions for use prior to use.
- This product is intended to be used by only medical specialist of fertility treatment.
- Aseptic technique should be used.
- Use only sterilized equipment and materials.
- Decontaminate the workroom.
- Follow procedures in an environmentally controlled room.
- The long term safety is unknown.

Copyright © KITAZATO CORPORATION All Rights Reserved.

- ### References
- Cobo A. Obstetric and perinatal outcome of babies born from vitrified oocytes. Fertility and Sterility, 2014.
 - Renzi L., Consistent and predictable delivery rates after oocyte vitrification: an observational longitudinal cohort multicentric study. Human Reproduction, 2012.
 - Shi W., Perinatal and neonatal outcomes of 494 babies delivered from 972 vitrified embryo transfers. Fertility and Sterility, 2012.
 - Cobo A., Outcomes of vitrified early cleavage-stage and blastocyst-stage embryos in a cryopreservation program: evaluation of 3,150 warming cycles. Fertility and Sterility, 2012.
 - Trokudes KM., Comparison outcome of fresh and vitrified donor oocytes in an egg-sharing donation program. Fertility & Sterility, 2011.
 - Inoue F., Hydroxypropyl cellulose as a macromolecular supplement for cryopreservation by vitrification of bovine oocytes and blastocysts and human oocytes.ESHRE and ASRM, 2011.