# Instructions For Use

# Cryotop®SC

Ver.6

## EN : INTENDED USE

CryotopSC is used for storage of vitrified human oocytes and embryos.

## QUALITY CONTROL TESTING

Each lot of CryotopSC receives the following tests:

- Je par la current USP Sterility Test
- Endotoxin by LAL methodology
- Mouse Embryo Assay (one cell)

All results of each lot are reported on a Certificate of Analysis which is available upon request.

## STORAGE INSTRUCTIONS AND STABILITY

Store in original sterile pack at 15-30 degree Celsius.

CryotopSC is stable until the expiration date shown on package label when stored as directed.

## PRECAUTIONS AND WARNING

Sterile if the package is unopened or undamaged. Do not use if package is broken.

Each CryotopSC is intended for single use.

Perform all vitrification and thawing operations for oocytes or embryos at room temperature (23-27 degree Celsius).

Liquid nitrogen may cause freeze burns if in contact with skin or eyes. Use appropriate precautions when work with liquid nitrogen.

Do not re-use. Re-use may cause change in product quality and raise the risk of poor embryo survival rate.

This product is intended to be used by medical specialists trained in fertility treatment. Aseptic technique should be used.

## INSTRUCTION FOR USE

### Vitrification

1. Open the sterile pack of CryotopSC and write necessary information about a patient on the handle of CryotopSC.
2. Fill 90% of a container with fresh liquid nitrogen.
3. Prepare oocytes or embryos for vitrification according to laboratory protocol.
4. Place the CryotopSC under the microscope and adjust the focus on the black mark of the CryotopSC. Logo should be up.
5. Gently place the oocytes or embryos by the black mark of CryotopSC sheet with minimal volume of vitrification solution using a suitable pipette.

**NOTE:** The recommended load of the CryotopSC is up to 4 oocytes or embryos.

6. Quickly plunge the CryotopSC tip to the middle of handle into liquid nitrogen. Then prop it against the edge of container.

7. Plunge the straw cap into liquid nitrogen and cut the above marking point.

8. Hold up the CryotopSC sheet within the 2.5cm height from the surface of liquid nitrogen. Insert the CryotopSC into the cap straw and seal it with a heat sealer.

**NOTE:** Do not let liquid nitrogen into the cap straw.

9. Immerse the capped CryotopSC into liquid nitrogen.

10. Put the CryotopSC in a cane and store it in a tank for long term storage.

**CAUTION:** Take care that the CryotopSC remains immersed in liquid nitrogen at all times until thawing.

### Thawing

1. Prepare the thawing media according to laboratory protocol.
2. Retrieve the cane which has the specific CryotopSC and quickly immerse the cane in a container filled with fresh liquid nitrogen. Collect the specific CryotopSC from the cane in the liquid nitrogen. Check the information about the patient on the handle of the CryotopSC.

**CAUTION:** The CryotopSC except the handle remain immersed in liquid nitrogen at all time. Place a container for liquid nitrogen by the stereo microscope.

3. Cut the sealed part of the straw cap with scissors remaining the tip of CryotopSC in liquid nitrogen.
4. Remove the CryotopSC from the straw cap, and quickly immerse the CryotopSC sheet into the thawing solution. It should be one second.
5. One minute after immersing into the thawing solution, gently aspirate the oocytes or embryos using a suitable pipette.
6. Perform the thawing procedure according to laboratory protocol.

## SYMBOLS

	Sterilized using irradiation		Use by
	Do not reuse		Storage Temperature
	Catalog Number		CE mark
	Batch Code		Do not use if package is damaged
	Consult Instructions for use		Do not re-sterilize
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community		

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