VORWORT

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THE IMPACT OF THE SHORT-TERM HUMAN SPERM STORAGE IN THE CRYOPROTECTANT-FREE MEDIUM ON SPERM MOTILITY AND VITALITY.

Nabil Sayme, Dr. med.¹, Marija Kijacic, Master of Biology Science², Thomas Krebs, Biology¹ and Dieter Maas, Prof. Dr. med.¹,
(1)Team Kinderwunsch Hannover, Hannover, Germany, (2)Saarland University Medical Center, Homburg, Germany

ABSTRACT

Slow freezing is currently the most commonly used technique for sperm cryopreservation since the vitrification of spermatozoa is still a rather unexplored methodology. Storage sperm at +4 °C is a relatively new technique and previously shown that cool storage without cryoprotectants is better for DNA integrity. In the aim of reach better recovery rates, many studies confirmed that freezing/store sperm without cryoprotectants gives better results.

OBJECTIVES

The purpose of the study was to investigate does it the sperm storage in the cryoprotectant-free medium a good alternative for short-time preservation (up to 4 weeks) compare to the conventional slow freezing, as well as the impact of the short-term human sperm storage on sperm motility and vitality. The study included 20 sperm samples collected between February 2018 - April 2018. Out of 20 samples, 10 were normozoospermic, 5 were teratozoospermic and 5 were asthenozoospermic.

METHODS

Native samples with volume higher than 3 ml before preparation was divided equally in the aim to reach the same sperm concentration of 20 Mill/ml. After gradient preparation, the volume of 200µl sample was treated with the same volume of cryoprotectant (GM501 Sperm store, Gynemed, Germany) or sperm preserve medium (Sperm Preserve, Gynemed, Germany). After slow freezing two straws of each sample were preserved into liquid nitrogen for a period of one and four weeks. Samples treated with Sperm Preserve were divided as well into two straws and stored in fridge on 4 °C degrees for the same period. After this period samples were thawed with Sperm Active (GM501 Sperm Active, Gynemed, Germany) and motility and vitality after these two freezing procedures were compared. For statistical analysis, a One-Way ANOVA was used.



Fig. 2

RESULTS

After one week of slow freezing and storage in liquid nitrogen, 24.9±11.22% of spermatozoa regained their motility compare to samples which were stored on 4 °C where recovery rate was 32.75±11.02% (Fig. 1A). One-Way ANOVA confirmed that there is a significant difference between these two groups (p= 0.31). The sperm motility rate after four weeks was slightly lower 20.4±6.87% in the slow freezing sample group, compare to 28.05±9.81% after storage at 4 °C but still, further statistical analysis confirmed a significant difference between these two groups (p= 0.005). Asthenozoospermic samples stored at 4 °C had better motility recovery rate after one week 28.8±11.6% vs 13.6±2.96% (p=0.02) than after four weeks 24.2±14% vs 12.8 ±4.7% where the difference between these two groups was not statistically significant as well as neither between teratozoospermic samples. Vitality was one of the characteristics which we analyzed (Fig. 1B) as well and the difference was significant especially after one week (p= 0.0001) where survival rate after slow freezing was 40.5±11.80% compared to the storage sample where that number was 54.25±13.20 %. After four weeks as well, a higher percentage of sperm survive in the storage group 39.75±7.88% compare to the 30.7±6.78% of slow frozen samples (p=0.001).

CONCLUSIONS

The cryoprotectant-free sperm storage protocol tested in this study renders considerably better recovery rates (motility and vitality) of the sperm compared to slow freezing and could be an option for short-term storage of human semen

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Herausgeber: GYNEMED GmbH & Co. KG