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Intended use

Injection pipettes and holding pipettes

- The products are designed for usage in the field of Assisted Reproductive Techniques (ART) for Intracytoplasmic Sperm injection (ICSI). With the holding pipette an oocyte is held by suction. With the injection pipette a spermatozoon is aspirated and injected into the oocyte using a micromanipulator under microscopic control.

Blastomer biopsy pipettes

- The blastomer biopsy pipette is used for the biopsy of blastomeres, as soon as the embryo is at the stage of 6-8 cells, or trophectoderm cells (TE). Blastomere and TE-cells are used for Preimplantation Diagnostics (PID). The utilization is done using a micromanipulator under microscopic control.

Product description

- All pipettes are made of borosilicate glass. For sterilization the products are gamma-irradiated and single-packed.
- The holding pipettes are available in one version of the tip (blunt). They can be ordered with different outer diameters (OD) and different bending angles (BA) at the tip.
- The injection pipettes are available with tips in versions with and without spike (spiked/bevelled). They can be ordered with different inner diameters (ID) and different bending angles (BA) at the tip.
- The blastomere biopsy pipettes are available with different versions of the tip (blunt/bevelled). They can be ordered with different inner diameters (ID) and different bending angles (BA) at the tip.
- The selection of the version of a pipette depends on the preferences and experiences of the user.

Pre-use checks


- Please check the package for damages. If the sterile package is damaged please dispose the product.

Preparation

- Grab the pipette between the two brackets of the plastic holder. As the first step release the pipette from anterior clip (tip) – otherwise the risk is involved to break the tip. Secondly release the pipette from the rear clip.
- Insert the micropipette in the micromanipulation holder and attach it corresponding to the instruction for use of your micromanipulator.

Instructions for use

- Hold the oocyte with the holding pipette.
- Aspirate a spermatozoon with the injection pipette and inject it into the oocyte using a micromanipulator.
- The blastomere biopsy pipette is used for the biopsy of blastomeres.

Instruction for Use Micropipettes

Quality control

- Products are sterilized by gamma-irradiation to a sterility assurance level (SAL) of 10^{-6} .
- Not toxic to embryos - proven by Mouse Embryo Assay (MEA) with the following requirements: day 6 $\geq$ 80% blastocysts. Certificate of analysis available for each batch on request.

Contraindications

- There are no known contraindications regarding the usage of these products.


Storage

- We recommend to store the goods under clean and dry conditions, protected from sunlight.

Disposal

- Take the micropipette out of the micromanipulation holder of your micromanipulator.
- Dispose the micropipette and the package in an appropriate container in accordance with your national regulations.

Explanation of symbols

Precautions and warnings

- The pipettes should only be used by qualified personnel. The intended users are IVF professionals (lab technicians, embryologists or medical doctors).
- Handle the products with care. The pipettes are made of thin borosilicate glass and are fragile – they carry the risk of cuts and punctures.
- Never point pipettes towards persons – high pressure can lead to a blow-pipe-effect which may cause incorrectly mounted pipettes to fly out of the micromanipulation holder.
- Do not rinse with fluids – rinsing involves the risk of contamination.
- Do not reuse – the products are intended for single use only; reusing involves a high risk of cross-contamination.
- Do not resterilize – resterilization involves a high risk to impair the functionality.
- Do not use if the product or the package are damaged prior to usage – this would involve the risk of e.g. contamination.
- Protect from direct sunlight – UV-light compromises the functionality of the sterile system and can reduce the shelf life.
- Protect from humidity – humidity compromises the functionality of the sterile system and can reduce the shelf life.
- Do not use after expiry date – safety and functionality are proven up to that date but not after.
- Only made for intended use.