Material included in the package:

- Product code: 4 GM501 SpermMobil-5-S
- 1 x 1 ml GM501 SpermMobil
- Product code: 4 GM501 SMOBILS
- 1 x 5 ml GM501 SpermMobil

Material not included in the package:

- Incubator (no CO₂)
 - Petri-dishes
 - LAF bench (ISO 5 environment)
 - Microscope
 - Test tubes
- Product specifications and quality control:**
- All components are of highest available purity (European Pharmacopoeia and/or USP standard) if applicable.
 - A certificate of analysis is available for each batch upon request from our website with respective lot number.
 - The MSDS for GM501 SpermMobil is available upon request and can also be downloaded from our website.
 - GM501 SpermMobil is manufactured and tested according to the following specifications:

pH (at 37°C)	7.20-7.60
Osmolality (mOsm/kg)	280-320
Sterility	sterile - SAL10 ⁻³ (Sterility Assurance Level)
Endotoxins (EU/ml)	< 0.50

Sperm Survival test
≥80% survival after 4 hours exposure of density selected spermatozoa to the test medium in %

EN: Intended use/Intended users:

- GM501 SpermMobil is a HEPEs-buffered reagent containing low bi-carbonate. It is used for in-vitro examination of sperms of necrozoospermic ejaculates as well as of immotile sperms isolated from testicular tissue (TESE).
 - The intended users are IVF professionals (lab technicians, embryologists or medical doctors).
- Composition:**
- NaCl, KCl, KH_{PO}, MgSO₄, 7H₂O, NaHCO₃, CaCl₂H₂O, D(+)-Glucosaminohydrat, Na-laktat, Na-pivrat, EDTA, Alanyl-Glutammina, water, non-essential and essential Amino Acids, HEPEs, Theophylline, Phenol Red.
- Instructions for use:**
- Do not equilibrate GM501 SpermMobil in a CO₂-incubator, just warm it up to 37°C. GM501 SpermMobil is HEPEs-buffered. Incubation in a CO₂-incubator will lower the pH.
 - To facilitate spermactivation add 3-4 µl GM501 SpermMobil to the sperm suspension (30-40µl/ultracentrifugation 1:10) in sperm preparation-medium inside the Petri-dishes.
 - GM501 SpermMobil should be added to the opposing side from where the sperms are to be aspirated.
 - Wait for 10 minutes. The activating effect initiates after a few minutes and lasts for up to one hour.
 - The dish should be placed on a heating plate (37°C) during the diagnostic evaluation.
- Precautions and warnings:**
- Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.
 - Always work under strict hygiene conditions (ISO 5 environment, e.g. LAF-bench) to avoid potential contamination.
 - Only for the intended use.
 - The user of this device is responsible for maintaining traceability of the product and must comply with national regulations regarding traceability where applicable.
- Pre-use instructions:**
- Do not use the product if it becomes discoloured, cloudy, or shows any evidence of microbial contamination.
 - Do not use the product if seal of the container is opened or defective when the product is delivered.
 - The concentration of the xanthine derivative Theophylline leads to the occasional formation of precipitate at 2-8°C storage temperature. It is dissolved completely when warmed at 37°C (within 24 h) and does not compromise function and safety of the product.

Storage instructions and stability:

- The shelf life is 6 month from time of manufacture.
- Store between 2-8°C.
- Do not freeze before use.
- Keep away from sun (light).
- The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C). Content cannot be re-sterilized after opening.
- Do not use after expiry date.

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