

# Nature's Principles. Proven Success.


## Understanding ZyMöt® Sperm Separation Devices

### Better Sperm Selection by Mimicking Nature


ZyMöt Fertility, Inc., has developed novel devices for use in ART procedures conducted by fertility clinics and OB/GYN practices. ZyMöt Multi Sperm Separation Devices prepare motile sperm from semen for use in IUI, IVF and ICSI procedures. These devices are the first of their kind and available worldwide.

Our revolutionary tools simulate the natural barriers of the cervical and uterine pathway that sperm must overcome to fertilize an egg. We enable separation of optimally functional sperm without the use of damaging chemicals or density gradient centrifugation (DGC).

### ZyMöt Multi Device



The ZyMöt Multi device is available in two processing volumes, 850µL and 3mL. A sample is applied through the device's inlet port, connected to a lower sample chamber. This chamber is separated from an upper collection chamber by an 8µm microporous filter. Filter size was determined after comparison between 3µm, 5µm and 8µm pore sizes. Incubation times of 15, 30 and 45 minutes were evaluated, with sperm saturation achieved at 30 minutes.<sup>1</sup> These parameters yielded optimal sperm collection efficiency and motility, with the 8µm pore demonstrating the highest degree of normal morphology.<sup>2</sup>



During sample incubation, the most motile sperm migrate upward through the filter, leaving less motile sperm behind. Separated sperm are then collected from the upper chamber for subsequent use in IUI, IVF and ICSI procedures.

### Simplifying and Standardizing Workflow

Easy to adopt and simple to use, ZyMöt Multi Sperm Separation Devices provide considerable time savings and standardization over traditional methods. ZyMöt devices avoid damaging DGC and preserve normal sperm morphology. Contact us for more information about how to evaluate ZyMöt devices in your clinic. We offer comprehensive support with experts who are ready to help you incorporate our tools into your practice and extend your success. [Learn more at zymotfertility.com](https://www.zymotfertility.com).

### References

1. Tasoglu, S., Safaee, H., Zhang, X., Kingsley, J. L., Catalano, P. N., Gurkan, U. A., Nureddin, A., Kayaalp, E., Anchan, R. M., Maas, R. L., Tüzel, E. and Demirci, U. (2013), Exhaustion of Racing Sperm in Nature - Mimicking Microfluidic Channels During Sorting. *Small*, 9: 3374-3384. doi:[10.1002/smll.201300020](https://doi.org/10.1002/smll.201300020)
2. Asghar, W., Velasco, V., Kingsley, J. L., Shoukat, M. S., Shafiee, H., Anchan, R. M., Mutter, G. L., Tüzel, E. and Demirci, U. (2014), Selection of Functional Human Sperm with Higher DNA Integrity and Fewer Reactive Oxygen Species. *Adv. Healthcare Mater.*, 3: 1671-1679. doi:[10.1002/adhm.201400058](https://doi.org/10.1002/adhm.201400058)

# Revolutionizing Sperm Preparation

ZyMöt® Sperm Separation Devices: Better for your patients. Better for your practice.

## A Better Way to Prepare Sperm

Quality, accuracy and efficiency are central to the success of a fertility practice. Traditional sperm preparation methods are not only time-consuming and laborious, but could cause additional sperm DNA fragmentation<sup>1</sup> and cellular stress,<sup>2</sup> lowering the odds of success.<sup>3,4</sup> ZyMöt Sperm Separation Devices are a better way to prepare sperm for use in IUI, IVF and ICSI procedures. It's that simple.

## Simple to Adopt. Easy to Use.


Available worldwide, ZyMöt devices efficiently isolate healthy, rapidly-progressive sperm, to help achieve outcomes that matter.<sup>5,6</sup> Minimal training is required, with simple, standardized procedures that help users **quickly achieve optimal performance**.

## Work on Your Timeline

ZyMöt devices enable processing whenever a sample is ready, eliminating delays caused by an equipment bottleneck. With **only 5 minutes of total hands-on tech time per sample**, every ZyMöt-processed specimen represents a significant time savings over traditional methods. Using ZyMöt devices frees staff for other critical tasks and improves lab productivity.

## Fewer Steps. More Confidence.

A shorter chain of custody – fewer movements per sample – means that ZyMöt devices help minimize mismatching risk, reducing the potential for costly error. **Processing sperm with ZyMöt devices gives providers more confidence and gives patients more peace of mind.**



Comparison of major sperm separation steps when using ZyMöt Sperm Separation Devices (top) versus using the traditional method (bottom). ZyMöt requires fewer movements per sample, improving efficiency and productivity while reducing risk of costly errors.

## Natural. Simple. Effective.

Try ZyMöt Sperm Separation Devices and realize immediate savings of time and resources, while providing premium quality sperm separation for your patients. **Learn more at [zymotfertility.com](http://zymotfertility.com).**

## References

1. Zini, A., *et al.* Urology (2000). doi: 10.1016/S0090-4295(00)00770-6
2. Aitken, R. J. and Clarkson, J. S. J Andrology (1988). doi: 10.1002/j.1939-4640.1988.tb01067.x
3. Dimakopoulou A., *et al.* J Endocr Soc. (2019). doi: 10.1210/je.2019-OR18-5
4. Jayasena C.N., *et al.* Clin Chem (2019). doi: 10.1373/clinchem.2018.289348
5. Parrella, A., *et al.* J Assist Reprod Genet (2019). doi: 10.1007/s10815-019-01543-5
6. Gode, F. *et al.* Fertil Steril (2019). doi: 10.1016/j.fertnstert.2019.06.037

# Publication Spotlight: Utilizing Spermatozoa with Higher Genomic Competence Improves ICSI Outcomes

Understanding the latest science in the ZyMöt® revolution

## ZyMöt Device Helps Select Sperm with High Genomic Integrity

Building on research from Parrella *et al* in 2019<sup>1</sup>, scientists at Weill-Cornell Medical College have been examining various impacts of genomic integrity. Specifically, these Weill-Cornell researchers have studied sperm chromatin fragmentation (SCF) and the effects on clinical outcomes.

In a study presented at ASRM 2021, Keating's objective was to demonstrate that selecting spermatozoa with the highest genomic integrity utilizing the ZyMöt Multi 850µl Sperm Separation Device would enhance ICSI outcomes. In this follow-on to the Parrella *et al* 2019 study<sup>1</sup>, the first step was to look

**FIG. 1 SEMEN PARAMETERS (N=126 PATIENTS)**

	SELECTION METHOD		
	RAW	DGC	ZYMÖT
MOTILITY (%)	33.7±14	61.2±33	96.3±13
MORPHOLOGY (%)	2.2±1	2.0±1	3.0±1
SCF (%)	23	18	1.4

at sperm parameters in men known to have SCF >15% after processing by, in parallel, density gradient centrifugation (DGC) and the ZyMöt device.

The SCF results collected from these studies, as shown in **Fig. 1**, highlighted the need to continue to focus on genomic integrity for male-factor patients. In particular, the encouraging results using the ZyMöt Multi 850 device supported taking a further look at clinical outcomes in additional patients.


## Low SCF with ZyMöt Device Leads to Better Clinical Outcomes vs. DGC

Based on the significant improvement in SCF utilizing the ZyMöt Multi 850 device, the researchers next looked at the clinical outcomes of this device compared to DGC (see **Figs. 2-3**). A total of 21 men (aged 43.3±8 years) had an average SCF in their raw semen of 22.1±10%, which decreased to 19.1±7% after DGC sperm preparation. These men underwent 39 ICSI cycles with their female partners (aged 38.0±4 years).

**FIG. 2 SEMEN PARAMETERS (N=21 PATIENTS)**

	SELECTION METHOD		
	RAW	DGC	ZYMÖT
SCF (%)	22.1±10	19.1±7	1.2±1

**FIG. 3 CLINICAL OUTCOMES - FRESH EMBRYO TRANSFER**




Subsequently, these couples underwent 26 ICSI cycles utilizing sperm preparation with the ZyMöt Multi 850 device. The SCF after ZyMöt device use was 1.2±1%, substantially lower than the raw sample and DGC. All cycles were fresh embryo transfers<sup>2</sup>.

## ICSI Treatment with ZyMöt Device for Men with High SCF Enhances Clinical Outcomes

After seeing the significant results comparing the ZyMöt device to DGC, Keating sought to evaluate the clinical outcomes solely using the ZyMöt device in couples where the men were known to have elevated SCF. Fifty-five (55) men (aged 42.3±8 years) were treated in 69 ICSI cycles with their female partners (aged 38.3±5 years). The SCF in their raw samples was 22.3±10%, which fell to 3.0±4% (p < 0.0001) following ZyMöt sperm sample preparation (see **Fig. 4**).

As Keating *et al* concluded, "Compared to the more conventional DGC, MFSS (the ZyMöt device) is capable of selecting the most progressively motile spermatozoa with the highest genomic integrity. Treatment by ICSI with MFSS (ZyMöt) for men with high sperm DNA fragmentation enhances fertilization, embryo development, and clinical pregnancies."

**FIG. 4 CLINICAL OUTCOMES - ZYMÖT DEVICE ONLY**



## References

1. Parrella, A., Keating, D., Cheung, S. *et al*. J Assist Reprod Genet 36, 2057-2066 (2019), doi: 10.1007/s10815-019-01543-5
2. Keating, D., Tavares, D., Rosenwaks, Z., Palermo, G. Fertility and Sterility Vol. 116, E67-E68(2021). doi: <https://doi.org/10.1016/j.fertnstert.2021.07.190>



# Publication Spotlight: Utilizing Spermatozoa with Higher Genomic Competence Improves ICSI Outcomes

Understanding the latest science in the ZyMöt® revolution

## ZyMöt Device Helps Select Sperm with High Genomic Integrity

FIG. 1 SEMEN PARAMETERS (N=126 PATIENTS)


	SELECTION METHOD		
	RAW	DGC	ZYMÖT
MOTILITY (%)	33.7±14	61.2±33	96.3±13
MORPHOLOGY (%)	2.2±1	2.0±1	3.0±1
SCF (%)	23	18	1.4

## Low SCF with ZyMöt Device Leads to Better Clinical Outcomes vs. DGC

FIG. 2 SEMEN PARAMETERS (N=21 PATIENTS)


	SELECTION METHOD		
	RAW	DGC	ZYMÖT
SCF (%)	22.1±10	19.1±7	1.2±1

FIG. 3 CLINICAL OUTCOMES - FRESH EMBRYO TRANSFER



## ICSI Treatment with ZyMöt Device for Men with High SCF Enhances Clinical Outcomes

FIG. 4 CLINICAL OUTCOMES - ZYMÖT DEVICE ONLY



# Publication Spotlight: A Sperm Selection Technique to Improve Embryo Ploidy

Understanding the latest science in the ZyMöt® revolution

## OBJECTIVE

To assess the role of an enhanced sperm selection method in mitigating paternal contributions to embryo aneuploidy.

## MATERIALS AND METHODS

Over the last 4 years, 57 couples underwent ICSI with sperm selected by density gradient centrifugation (DGC), resulting in few frozen embryo transfers (FETs) due to consistent embryo aneuploidy following preimplantation genetic testing for aneuploidy (PGT-A).

These men consented to sperm chromatin fragmentation (SCF) assessment, inclusive of double-stranded DNA breaks (dsDNA) in their raw semen, as well as post-DGC and post-microfluidic sperm selection with the ZyMöt Multi Sperm Separation Device (ZyMöt). These couples underwent subsequent ICSI cycles with ZyMöt. Outcomes of cycles processed by DGC and ZyMöt were analyzed and compared.

SCF was assessed by terminal deoxynucleotidyl dUTP transferase nick-end labeling (TUNEL) on  $\geq 500$  spermatozoa per patient, with a normal threshold of  $\leq 15\%$ . A neutral Comet assay was used to assess dsDNA on  $\geq 200$  spermatozoa, utilizing a modified in-house protocol and a normal threshold of  $\leq 3\%$ .

## RESULTS

A total of 57 men had the following semen parameters: concentration of  $40.0 \pm 32 \times 10^6/\text{mL}$ ,  $37.1 \pm 11\%$  motility, and  $2.2 \pm 1\%$  normal morphology. After selection by DGC or ZyMöt, the concentrations were  $3.3 \pm 3.4 \times 10^6/\text{mL}$  and  $8.0 \pm 13 \times 10^6/\text{mL}$ , with  $58.0 \pm 29\%$  and  $96.9 \pm 9\%$  motility, respectively ( $p < 0.0001$ ). The SCF decreased from  $21 \pm 14\%$  in raw specimens to  $18 \pm 6\%$  following DGC and to  $1.9 \pm 1\%$  following ZyMöt ( $p < 0.001$ ). The dsDNA fell from  $3.6 \pm 2\%$  in raw specimens to  $3.1 \pm 1\%$  after DGC and to  $0.3 \pm 0.2\%$  after ZyMöt ( $p < 0.001$ ).

These men (aged  $40.9 \pm 6$  years) underwent DGC selection for 71 ICSI cycles with their female partners (aged  $36.5 \pm 5$  years), achieving a fertilization rate of  $58.4\%$  ( $403/690$ ) and a blastocyst euploidy rate of  $28.5\%$  ( $47/165$ ). Only 19 FET cycles were performed, with a  $6.7\%$  implantation rate ( $2/30$ ) and two clinical pregnancies resulting in miscarriage.

Subsequently, these men had their specimens selected by ZyMöt in 71 ICSI cycles, resulting in a higher fertilization rate of  $75.9\%$  ( $647/852$ ;  $p < 0.0001$ ) and a much improved ( $p < 0.0001$ ) blastocyst euploidy rate of  $48.9\%$  by PGT-A ( $192/389$ ). In 48 FET cycles, 51 embryos were replaced with an increased implantation rate of  $60.8\%$  ( $31/51$ ;  $p < 0.0001$ ), a CPR of  $64.6\%$  ( $31/48$ ;  $p < 0.001$ ), and an ongoing/delivery rate of  $62.5\%$  ( $30/48$ ;  $p < 0.0001$ ).

FIG. 1 PATIENT DEMOGRAPHICS PGA-T (N=57 COUPLES)

CYCLES	SELECTION METHOD	
	DENSITY GRADIENT	ZYMÖT DEVICE
MATERNAL AGE (M $\pm$ SD)	71	71
MATERNAL AGE (M $\pm$ SD)	36.5 $\pm$ 5	37.2 $\pm$ 15
MATERNAL AGE (M $\pm$ SD)	40.9 $\pm$ 6	41.1 $\pm$ 7

FIG. 2 EUPLOIDY RATES - PGT-A





FIG. 3 CLINICAL OUTCOMES



## CONCLUSIONS

With its dsDNA component, SCF tangibly contributes to embryo structural chromosomal abnormalities. An enhanced spermatozoa selection method for ICSI appears to remarkably increase the proportion of euploid blastocysts with consequent successful clinical outcomes.

## IMPACT STATEMENT

Sperm genomic integrity is associated with the ploidy of the conceptus, and a high SCF inclusive of dsDNA can be mitigated by proper sperm selection.

## References

- Keating, D., Tavares, D., Rosenwaks, Z., Palermo, G. Fertility and Sterility Vol. 116, P-53(2021). doi: <https://doi.org/10.1016/j.fertnstert.2021.07.380>

# Publication Spotlight: DNA & ROS Levels After ZyMöt® Sperm Prep

Understanding the latest science in the ZyMöt revolution

## The Need for Healthy Sperm


Using healthy sperm for IUI, IVF and ICSI procedures is more important than ever. In new research<sup>1</sup> from scientists at Imperial College London, recurrent pregnancy loss was directly connected to the presence of elevated sperm DNA fragmentation and reactive oxygen species, along with a lower percentage of normal morphology. This follows a growing body of evidence that links improved sperm health to better pregnancy outcomes.<sup>2</sup> ZyMöt Sperm Separation Devices enable the preparation of sperm with low levels of DNA fragmentation and oxidative stress.

## Avoiding DNA Fragmentation and Oxidative Stress

ZyMöt devices have been shown to separate sperm with near-zero DNA fragmentation, compared to density gradient centrifugation.<sup>3</sup> In an independent study from Midwest Fertility Specialists, ZyMöt Multi (850µL) Sperm Separation Devices were directly compared to traditional sperm preparation techniques.<sup>4</sup> This clinical research determined which approach resulted in improved DNA fragmentation index (DFI) and other sperm health biomarkers such as oxidative stress adducts (OSA) and high DNA stainability (HDS).

**Results:** Using ZyMöt devices significantly reduced DFI ( $P < 0.05$ ) compared to standard protocols: two commercially available gradients, and gradients followed by swim-up. The device also effectively reduced ( $P < 0.05$ ) OSA levels, a measurement of oxidative stress, and HDS, a measurement of immature cells and high histone retention. "Overall, the quality of the sperm obtained post-processing was improved by the use of the separation device," wrote the study author.

**Conclusion:** Using ZyMöt devices shows statistically significant improvements in three DNA- and stress-focused indicators of sperm health and function, when compared to traditional methods.



## Improving Efficiency and Outcomes

ZyMöt devices are easy to adopt and simple to use, helping labs quickly achieve optimal performance. With only 5 minutes of total hands-on tech time per sample, every ZyMöt-processed specimen represents a significant time savings over traditional methods. In addition to increased efficiency, ZyMöt devices deliver improved sperm performance to help achieve the best possible outcomes in IUI, IVF and ICSI procedures. Learn more at [zymotfertility.com](https://zymotfertility.com).

## References

1. A.P. Dimakopoulou *et al.* Elevated semen oxidative stress in male partners as novel marker of recurrent pregnancy loss. ENDO 2019, New Orleans, March 24, 2019.
2. C.N. Jayasena *et al.* Reduced testicular steroidogenesis and increased semen oxidative stress in male partners as novel markers of recurrent miscarriage. Clinical Chemistry. Volume 65, January 2019, p. 161.
3. M.M. Quinn *et al.* Microfluidic sorting selects sperm for clinical use with reduced DNA damage compared to density gradient centrifugation with swim-up in split semen samples. Hum Reprod. July 10, 2018. doi: 10.1093/humrep/dey239.
4. Broussard A *et al.* Sperm DNA fragmentation (SDF) was most effectively improved by a sperm separation device compared to different gradient and swim-up methods. Fertility and Sterility, Volume 111, Issue 4, e15.






# Publication Spotlight: Improving Genomic Integrity & Patient Outcomes

Understanding the latest science in the ZyMöt® revolution

## Not All Sperm are Equal

Using the best sperm helps increase the odds of a successful fertility treatment cycle. Not all sperm are created equal: up to 11% of men with a “normal” semen analysis have a measurable problem with sperm chromatin (DNA) fragmentation, and thus reduced motility.<sup>1</sup> Double-stranded sperm DNA damage is a cause of delay in embryo development and can impair implantation rates.<sup>2</sup>

ZyMöt Sperm Separation Devices are a new way to process sperm. ZyMöt devices enable the separation of sperm with nearly undetectable levels of DNA fragmentation and oxidative stress.<sup>3</sup> Improved sperm health means better clinical outcomes.<sup>4-6</sup>



## Results: Sperm DNA Damage Lowers the Odds of Success

Work from Keating showed that double-stranded DNA breaks in sperm were a major factor in chromosomal abnormalities, embryo aneuploidy and pregnancy loss.<sup>7</sup> This highlights the need to focus on genomic integrity – not just for male-factor patients, but for every sample.

## Results: Improved Outcomes for Challenging Patients

In an update to her 2019 publication<sup>1</sup>, Parrella and colleagues studied patients with histories of ART failure and high sperm chromatin fragmentation (SCF) ( $\geq 22\%$ ). This research asked if microfluidic sperm separation was able to select sperm with higher chromatin integrity.<sup>8</sup>

One patient group underwent fresh embryo transfer (FET) after processing with DGC. Initially, this group saw low levels of clinical pregnancy and high levels of loss. These patients then had their semen specimens processed with ZyMöt in a subsequent ICSI cycle, yielding significantly higher implantation rates, clinical/ongoing pregnancy rates, and decreased pregnancy loss.

In another group, patients also had both high SCF and a history of high embryo aneuploidy rates. Patients underwent PGT-A and frozen embryo transfer, after sperm processing with either DGC or the ZyMöt device. Euploidy rates were significantly higher with the ZyMöt device compared to DGC processing. Implantation rates, clinical pregnancy rates, ongoing/delivered rates (there were none with DGC) were all significantly higher with ZyMöt compared to DGC processing with greater pregnancy loss, respectively.

## The ZyMöt Difference

The science is clear: it's essential to do everything we can to improve sample quality by selecting sperm with the lowest possible levels of DNA fragmentation. Processing with the ZyMöt device enhances sperm sample motility, and progression and morphology, along with providing a “remarkable reduction” of DNA fragmentation.<sup>1</sup> ZyMöt devices yield sperm with higher genomic competence, demonstrated by their improved euploid rate and ability to establish healthy pregnancies, even for couples with histories of previous ART failure. Learn more at [zymotfertility.com](https://zymotfertility.com).

## References

1. Parrella, A. *et al.* *J Assist Reprod Genet* (2019). doi: 10.1007/s10815-019-01543-5
2. Cassanova, A. *et al.* *Fertility and Sterility* (2019). doi: 10.1016/j.fertnstert.2018.11.035
3. Broussard, A., *et al.* *Fertility and Sterility* (2019), Volume 111, Issue 4, e15
4. Dimakopoulou A., *et al.* *J Endocr Soc.* (2019). doi: 10.1210/js.2019-OR18-5
5. Jayasena C.N., *et al.* *Clin Chem* (2019). doi: 10.1373/clinchem.2018.289348
6. Palermo, G. *et al.* *Fertility and Sterility* Volume (2014). doi: 10.1016/j.fertnstert.2014.10.008
7. Keating, D. *Fertility and Sterility* (2020), Volume 114, Issue 3, P-496
8. Parella A., *et al.* *Fertility and Sterility* (2020), Volume 114, Issue 3, O-69 & P-781

# Publications Spotlight: Comparison of Microfluid Sperm Sorting Chip and Density Gradient Methods for Use in Intrauterine Insemination Cycles

Understanding the latest science in the ZyMöt® revolution

## A Better Way to Process Sperm for IUI


ZyMöt Sperm Separation Devices offer an alternative preparation method that allows for simple, natural, and effective isolation of motile sperm.

**Results:** In a retrospective study of 265 IUI patients with unexplained infertility, patients whose semen samples were processed using ZyMöt Sperm Separation Devices were 3.5 times more likely to achieve an ongoing pregnancy than the age-matched control group, where sperm was processed with the traditional, centrifugation-based method.<sup>1</sup> In this study, the ZyMöt IUI treatment group also experienced a reduced miscarriage rate when compared to patients whose semen samples were processed by density gradient (0% vs. 5% respectively).<sup>1</sup>

**Conclusion:** Processing sperm with ZyMöt devices can improve a patient's IUI treatment prognosis at the onset of their infertility journey, offering a greater chance of success with a lower risk treatment option.



## ONGOING PREGNANCY RATE



The effect of IUI sperm preparation method on pregnancy outcomes: ZyMöt improved treatment prognosis.

## Simple to Adopt. Easy to Use.

FDA-cleared, CE-certified and available worldwide, ZyMöt Sperm Separation Devices efficiently isolate the healthiest, rapidly progressive sperm, to help achieve outcomes that matter.<sup>2</sup> Minimal training is required, with simple, standardized procedures that help users quickly achieve optimal performance. ZyMöt Sperm Separation Devices are a better way to prepare sperm. It's that simple. Learn more at [zymotfertility.com](https://zymotfertility.com).

## References

1. Gode, F., et al. Fertil Steril (2019). doi: 10.1016/j.fertnstert.2019.06.037
2. Parrella, A., et al. J Assist Reprod Genet (2019). doi: 10.1007/s10815-019-01543-5



# Publications Spotlight: Improved Patient Outcomes after ZyMöt® Device Prep

Understanding the latest science in the ZyMöt revolution

## Recent Studies Examined Euploidy and Ongoing Pregnancy Rates

ZyMöt Sperm Separation Devices have been designed and developed to aid reproductive medicine professionals in the selection of motile sperm for use in assisted reproductive technology (ART) procedures. In new research presented at ASRM 2020, investigators examined euploidy and ongoing pregnancy rates, and saw significant improvement when processing samples with ZyMöt devices.


## Results: Improved Euploidy Rates

Anderson and colleagues conducted a prospective cohort study<sup>1</sup> that compared the impact of sperm prepared utilizing density gradient centrifugation (DGC) or sperm separation with the ZyMöt Multi (850µL) device on euploidy and pregnancy outcomes. The D5 euploid rate was significantly higher using ZyMöt compared to DGC (below, left). Anderson also presented results based on a Six Sigma-style evaluation of time and showed that ZyMöt saves procedural steps and time.

Beyhan and colleagues conducted a retrospective study<sup>2</sup> that examined preimplantation development following ICSI after ZyMöt or DGC, in presumed normal to moderate male infertility patients. Similar fertilization and blastocyst conversion rates between the cohorts were observed. An increased euploid rate was observed for the ZyMöt-processed samples (below, middle).

## Results: Improved Ongoing Pregnancies

In another retrospective study, Palmerola and colleagues<sup>3</sup> compared ongoing pregnancy rates for two cohorts that used either DGC or ZyMöt device preparation. A significant improvement in ongoing pregnancies following single, euploid embryo transfer was observed (below, right). Fertilization, usable blastocysts and D5 and D6 biopsy rates were similar between the DGC and ZyMöt groups.



## Improving Efficiency and Outcomes

ZyMöt devices are simple to use, helping labs quickly achieve optimal performance. With only 5 minutes of total hands-on tech time per sample, every ZyMöt-processed specimen represents a significant time savings over traditional methods. Learn more at [zymotfertility.com](https://zymotfertility.com).

## References

1. Anderson T., *et al.*, Fertility and Sterility (2020), Volume 114, Issue 3, O-104
2. Right: Beyhan Z., *et al.*, Fertility and Sterility (2020), Volume 114, Issue 3, P-48
3. Palmerola K., *et al.*, Fertility and Sterility (2020), Volume 114, Issue 3, P-45

A Vital Innovation  
for Fertility Patients.



# A GAME-CHANGER FOR YOUR LAB.

ZyMöt® Sperm Separation Devices


MORE THAN  
**50%**  
FEWER STEPS

ZyMöt devices require less than half the sample handling steps compared to average sperm prep, greatly reducing the risk for error.




### MORE EMBRYO TRANSFERS

Increased euploidy rates mean more opportunities for embryo transfers.



### ALL-TIME HIGH DEMAND

Patient demand for ZyMöt devices is at an all-time high. Making them available in your lab will help attract more patients.



### CONSISTENCY

Ensure you're always getting the best sample every time, regardless of tech experience.

UP TO  
**80%**  
GREATER EFFICIENCY



Reduce hands-on prep time by up to 80% and free up your lab technician's time for other tasks.

**90%**  
LESS WASTE



Reduce the amount of sperm preparation media, tubes, pipettes, and other equipment by up to 90%.

**90%**  
CUSTOMER  
RETENTION




More than 90% of labs that have used ZyMöt devices continue to use them.

### SCALABILITY


ZyMöt devices streamline a typically laborious process, which saves you time and positions your lab for growth.




ZyMöt Fertility, Inc. | a DxNow, Inc., subsidiary  
401 Professional Drive, Suite 130, Gaithersburg, MD USA 20879-3429  
240.454.9893 | zymotfertility.com | info@zymotfertility.com




**1** Allow sample to liquify




**2** Draw 850µL of the sample.




**3** a) Achieve seal. b) Slowly inject sample.



**4** a) Draw 750µL of media. b) Prime outlet channel. c) Cover membrane surface.



**5** Incubate at 37°C for 30 minutes.




**6** Slowly aspirate a maximum of 500µL

### ICSI and IUI

Transfer the collected sample to an appropriate culture tube: a 4mL round bottom culture tube with a snap top or into the bottom of a 15mL conical tube.

Tubes using HEPES-buffered media may be held on the benchtop or tightly capped in an incubator.

Tubes using bicarbonate-buffered media should be stored in a CO<sub>2</sub> incubator with the lid loosely closed.



### IVF

Transfer the collected sample into a 15ml conical tube.

Add 3ml of bicarbonate-containing media (whatever media is usually used for the final suspension of sperm for conventional insemination) to the conical tube. Mix gently.

Centrifuge the conical tube for 5min at 300 x g.

Remove the supernatant, being careful to not disturb the lower pellet.



Perform count and motility as usual and dilute if needed to achieve appropriate final insemination concentration.

Store tube in a CO<sub>2</sub> incubator until insemination.

Insemination should occur more than 1hr, but less than 4hrs after preparation.



**7** Sample Handling After Collection – ICSI and IUI

**7** Sample Handling After Collection – IVF





**1** Allow sample to liquify

**2** Draw 3mL of the sample.



**3** a) Achieve seal. b) Slowly inject sample.

**4** a) Draw 2.5mL of media. b) Prime outlet channel. c) Cover membrane surface.



**5** Incubate at 37°C for 30 minutes.


**6** Slowly aspirate a maximum of 1mL.

#### ICSI and IUI

Transfer the collected sample to an appropriate culture tube: a 4mL round bottom culture tube with a snap top or into the bottom of a 15mL conical tube.

Tubes using HEPES-buffered media may be held on the benchtop or tightly capped in an incubator.

Tubes using bicarbonate-buffered media should be stored in a CO<sub>2</sub> incubator with the lid loosely closed.



#### IVF

Transfer the collected sample into a 15ml conical tube.

Add 3ml of bicarbonate-containing media (whatever media is usually used for the final suspension of sperm for conventional insemination) to the conical tube. Mix gently.

Centrifuge the conical tube for 5min at 300 x g.

Remove the supernatant, being careful to not disturb the lower pellet.

Perform count and motility as usual and dilute if needed to achieve appropriate final insemination concentration.

Store tube in a CO<sub>2</sub> incubator until insemination.

Insemination should occur more than 1hr, but less than 4hrs after preparation.

**7** Sample Handling After Collection – ICSI and IUI

**7** Sample Handling After Collection – IVF