

**Certificate of Analysis  
GM501 SpermAIR**

Batchnumber	AIR-1811
Production date	23 November 2018
Expiry date	31 May 2019
Total volume (mL)	25640
Batchnumber HSA	B3AKC00053

	<b>Results</b>	<b>Specifications</b>
pH	7,4	7,20-7,50
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 11 December 2018

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*