

**Certificate of Analysis
GM501 SpermAIR**

Batchnumber	AIR-2309
Production date	15 September 2023
Expiry date	31 March 2024
Batchnumber HSA	B03H003853
Total volume (mL)	26960

Parameter	Results	Specifications
pH	7,46	7,20-7,50
Osmolality	282	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	02 October 2023
Released by	Nicole Geiger
Signature	

REV 07_00 / 2023-09-15