


**Certificate of Analysis
GM501 SpermAIR**

Batchnumber AIR-2408
 Production date 02 September 2024
 Expiry date 31 March 2025
 Batchnumber HSA B03H100413
 Total volume (ml) 22040

Parameter	Results	Specifications
pH	7,34	7,20-7,50
Osmolality	279	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/mL
Concentration HSA	4,8	4,0-6,0 g/l
Sterility test (performed according to the current Ph. Eur. 2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	95	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 18 September 2024
 Released by Nicole Geiger
 Signature



REV 09_00 / 2024-08-16