

**Certificate of Analysis  
GM501 SpermAIR**

Batchnumber AIR-2409  
 Production date 19 September 2024  
 Expiry date 31 March 2025  
 Batchnumber HSA B03H100413  
 Total volume (ml) 21280

Product was sterilized by aseptic processing techniques


Parameter	Results	Specifications
pH	7,34	7,20-7,50
Osmolality	282	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	4,7	4,0-6,0 g/l
Sterility test (performed according to the current Ph. Eur. 2.6.1 / USP <71>)	No growth	No growth
MEA (1-cell to expanded blastocyst within 96 hours)	90,48	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 09 October 2024

Released by Nicole Geiger

Signature



REV 09\_00 / 2024-08-16