

Certificate of Analysis GM508 CultActive

Batchnumber	CACT-1806
Production date	04 July 2018
Expiry date	31 January 2019
Total volume (mL)	990

	Results	Specifications
pH	7,38	7,00-7,60
Osmolality	421	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Release date 19 July 2018

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann