

**Certificate of Analysis  
GM508 CultActive**

Batchnumber	CACT-1905
Production date	16 May 2019
Expiry date	30 November 2019
Total volume (mL)	1000

	<b>Results</b>	<b>Specifications</b>
pH	7,32	7,00-7,60
Osmolality	417	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date 31 May 2019

Released by Dipl. Ing. Benthe Brauer

Signature

