

**Certificate of Analysis
GM508 CultActive**

Batchnumber	CACT-2003
Production date	13 March 2020
Expiry date	30 September 2020
Total volume (mL)	1000

	Results	Specifications
pH	7,48	7,00-7,60
Osmolality	422	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95	≥ 80 %

Release date 30 March 2020

Released by Vanessa Bierlich

Signature