

**Certificate of Analysis
GM508 CultActive**

Batchnumber	CACT-2102
Production date	19 February 2021
Expiry date	31 August 2021
Total volume (mL)	1500

	Results	Specifications
pH	7,33	7,00-7,60
Osmolality	424	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date	08 March 2021
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann