

Certificate of Analysis GM508 CultActive

Batchnumber	CACT-2103
Production date	24 March 2021
Expiry date	30 September 2021
Total volume (mL)	1500

	Results	Specifications
pH	7,4	7,00-7,60
Osmolality	426	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95	≥ 80 %

Release date	13 April 2021
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann