

Certificate of Analysis GM508 CultActive

Batchnumber	CACT-2105
Production date	21 May 2021
Expiry date	30 November 2021
Total volume (mL)	1500

	Results	Specifications
pH	7,49	7,00-7,60
Osmolality	416	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date	08 June 2021
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann