

Certificate of Analysis GM508 CultActive

Batchnumber	CACT-2110
Production date	01 October 2021
Expiry date	30 April 2022
Total volume (mL)	1500

	Results	Specifications
pH	7,2	7,00-7,60
Osmolality	420	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Release date	18 October 2021
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann