

**Certificate of Analysis
GM508 CultActive**

Batchnumber	CACT-2207
Production date	02 June 2022
Expiry date	31 December 2022
Total volume (mL)	2700

	Results	Specifications
pH	7,37	7,00-7,60
Osmolality	410	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date	20 June 2022
Released by	Vanessa Bierlich
Signature	

P.P. 