

**Certificate of Analysis
GM508 CultActive**

Batchnumber	CACT-2305
Production date	31 July 2023
Expiry date	31 January 2024
Total volume (mL)	2000

	Results	Specifications
pH	7,4	7,00-7,60
Osmolality	425	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Release date	16 August 2023
Released by	Nicole Geiger
Signature	