

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2306
 Production date 15 September 2023
 Expiry date 31 March 2024
 Total volume (ml) 2700

Parameter	Results	Specifications
pH	7,43	7,00-7,60
Osmolality	425	400-450 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/ml
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	100	≥ 80 %

Release date 02 October 2023
 Released by Nicole Geiger
 Signature 