

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2403
 Production date 13 March 2024
 Expiry date 30 September 2024
 Total volume (ml)

| Parameter | Results | Specifications |
|---|----------------|-----------------------|
| pH | 7,34 | 7,00-7,60 |
| Osmolality | 420 | 400-450 mOsm/kg |
| Bacterial screen | NO GROWTH | No growth |
| Fungal screen | NO GROWTH | No growth |
| Endotoxin | < 0,05 | < 0,25 EU/ml |
| Mouse Embryo Assay (1-cell to blastocyst within 96 hours) | 100 | ≥ 80 % |

Release date 28 March 2024
 Released by Nicole Geiger
 Signature

