

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2404
 Production date 22 April 2024
 Expiry date 31 October 2024
 Total volume (ml) 2000

Parameter	Results	Specifications
pH	7,39	7,00-7,60
Osmolality	427	400-450 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/ml
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	90	≥ 80 %

Release date 07 May 2024
 Released by Joana Sturm
 Signature

