

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2407
 Production date 01 August 2024
 Expiry date 28 February 2025
 Total volume (ml) 2000

Parameter	Results	Specifications
pH	7,35	7,00-7,60
Osmolality	419	400-450 mOsm/kg
Sterility test (performed according to the current Ph. Eu. 2.6.1 / USP <71>)	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/ml
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	86	≥ 80 %

Release date 19 August 2024
 Released by Joana Sturm
 Signature



REV 08_00 / 2024-08-16