

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2409
 Production date 01 October 2024
 Expiry date 30 April 2025
 Total volume (ml) 2700

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,33	7,00-7,60
Osmolality	422	400-450 mOsm/kg
Endotoxin	< 0,05	< 0,25 EU/ml
Sterility test (performed according to the current Ph. Eu. 2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	100	≥ 80 %

Release date 16 October 2024

Released by Nicole Geiger

Signature



REV 08_00 / 2024-08-16

