

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2501
 Production date 20 January 2025
 Expiry date 31 July 2025
 Total volume (ml) 2700

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,31	7,00-7,60
Osmolality	419	400-450 mOsm/kg
Endotoxin	< 0,05	< 0,25 EU/ml
Sterility test (performed according to the current Ph. Eu. 2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	100	≥ 80 %

Release date 04 February 2025

Released by Joana Sturm

Signature



