

Certificate of Analysis
GM508 CultActive

Batchnumber CACT-2503
Production date 01 April 2025
Expiry date 31 October 2025
Total volume (ml) 2700

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH (@37°C - 6% CO ₂)	7,29	7,00-7,60
Osmolality	419	400-450 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Sterility test (performed according to the current Ph. Eu. 2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	100	≥ 80 %

Release date 17 April 2025

Released by Nicole Geiger

Signature

pp. N. Geiger

REV 09_00 / 2025-03-19