

**Certificate of Analysis
GM508 CultActive**

Batchnumber CACT-2510
Production date 12 November 2025
Expiry date 31 May 2026
Total volume (ml) 2700

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH (@37°C - 6% CO ₂)	7,27	7,00-7,60
Osmolality	425	400-450 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Sterility test (performed according to the current Ph. Eu. 2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to blastocyst within 96 hours)	100	≥ 80 %

Release date 27 November 2025

Released by Joana Sturm

Signature

