

**Certificate of Analysis  
GM501 Cult with Gentamicin**

Batchnumber	CULT-1805G
Production date	21 June 2018
Expiry date	31 December 2018
Total volume (mL)	29000
Batchnumber HSA	M749A6671

	<b>Results</b>	<b>Specifications</b>
pH	7,23	7,20-7,45
Osmolality	284	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 09 July 2018

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*