

**Certificate of Analysis**  
**GM501 Cult with Gentamicin and Phenol Red**

Batchnumber	CULT-1806PRG
Production date	16 August 2018
Expiry date	28 February 2019
Total volume (mL)	6400
Batchnumber HSA	K809B6671

	<b>Results</b>	<b>Specifications</b>
pH	7,23	7,20-7,45
Osmolality	276	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	03 September 2018
Released by	Dr. Julia Heinzmann
Signature	

*p.p. Julia Heinzmann*