

## Certificate of Analysis

### GM501 Cult with Gentamicin and Phenol Red

Batchnumber	CULT-2005PRG
Production date	20 August 2020
Expiry date	28 February 2021
Batchnumber HSA	B3ALD01343
Total volume (mL)	5900

	<b>Results</b>	<b>Specifications</b>
pH	7,28	7,20-7,45
Osmolality	279	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date	04 September 2020
Released by	Dr. Julia Heinzmann
Signature	

*p.p. Julia Heinzmann*