

**Certificate of Analysis  
GM501 Cult with Gentamicin**

Batchnumber	CULT-2104G
Production date	28 May 2021
Expiry date	30 November 2021
Batchnumber HSA	K041A6871
Total volume (mL)	15200

	<b>Results</b>	<b>Specifications</b>
pH	7,2	7,20-7,45
Osmolality	281	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

*Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.*

Release date 14 June 2021

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*