

Certificate of Analysis

GM501 Cult with Gentamicin

Batchnumber	CULT-2111G
Production date	18 November 2021
Expiry date	31 May 2022
Batchnumber HSA	M106A6671
Total volume (mL)	22100

	Results	Specifications
pH	7,25	7,20-7,45
Osmolality	281	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	90,48	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	14 December 2021
Released by	Dr. Julia Heinzmann
Signature	

p.p. Julia Heinzmann