

**Certificate of Analysis
GM501 Cult with Gentamicin**

Batchnumber	CULT-2207G
Production date	23 September 2022
Expiry date	31 March 2023
Batchnumber HSA	K203D667G
Total volume (mL)	21100

	Results	Specifications
pH	7,31	7,20-7,45
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 18 October 2022

Released by Vanessa Bierlich

Signature

P.P. 