


**Certificate of Analysis
GM501 Cult with Gentamicin**

Batchnumber	CULT-2405G
Production date	10 July 2024
Expiry date	31 January 2025
Batchnumber HSA	K331A687M
Total volume (ml)	12900

Parameter	Results	Specifications
pH	7,25	7,20-7,45
Osmolality	280	270-290 mOsm/kg
Bacterial screen	No growth	No growth
Fungal screen	No growth	No growth
Endotoxin	< 0,05	< 0,25 EU/ml
HSA concentration	8,8	8,0-12,0 mg/ml
Mouse Embryo Assay ^(1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	01 August 2024
Released by	Nicole Geiger
Signature	



REV 08_00 / 2024-02-28