

Certificate of Analysis

GM501 Cult with Gentamicin and Phenol Red

Batchnumber	CULT-2406PRG
Production date	02 September 2024
Expiry date	31 March 2025
Batchnumber HSA	K331A687L
Total volume (ml)	6800
Product was sterilized by aseptic processing techniques	

Parameter	Results	Specifications
pH	7,29	7,20-7,45
Osmolality	280	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	10,1	8,0-12,0 g/l
Concentration gentamicin	11	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	90	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	18 September 2024
Released by	Nicole Geiger
Signature	

REV 09_00 / 2024-09-17