

**Certificate of Analysis
GM501 Cult with Gentamicin**

Batchnumber CULT-2604G
 Production date 23 March 2026
 Expiry date 30 September 2026
 Batchnumber HSA MP10A6871, KP28A6871
 Total volume (ml) 15200

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH (@37°C - 6% CO ₂)	7,26	7,20-7,45
Osmolality	281	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	9,7	8,0-12,0 g/l
Concentration gentamicin	13	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	90	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 07 April 2026
 Released by Dr. René Richarz
 Signature

