


Certificate of Analysis

GM501 Cult with Gentamicin and Phenol Red

Batchnumber	CULT-2610PRG
Production date	25 September 2026
Expiry date	31 March 2027
Batchnumber HSA	KP51A6872
Total volume (ml)	4800
Product was sterilized by aseptic processing techniques	

Parameter	Results	Specifications
pH (@37°C - 6% CO ₂)	7,31	7,20-7,45
Osmolality	283	270-290 mOsm/kg
Endotoxin level	< 0,05	< 0,25 EU/ml
Concentration HSA	10,1	8,0-12,0 g/l
Concentration gentamicin	10	5-15 mg/l
Sterility test (performed according to the current Ph.Eur.2.6.1 / USP <71>)	No growth	No growth
Mouse Embryo Assay (1-cell to expanded blastocyst within 96 hours)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date	13 October 2026
Released by	Nicole Geiger
Signature	

REV 10_00 / 2025-08-04