

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-1702
Production date	14 April 2017
Expiry date	31 October 2017
Total volume (mL)	32000

	Results	Specifications
pH	7,3	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date 02 May 2017

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann