

Certificate of Analysis GM501 Flush

| | |
|-------------------|-----------------|
| Batchnumber | FLUSH-1802 |
| Production date | 09 April 2018 |
| Expiry date | 31 October 2018 |
| Total volume (mL) | 21500 |

| | Results | Specifications |
|-----------------------------|---------|-----------------|
| pH | 7,26 | 7,20-7,40 |
| Osmolality | 278 | 270-290 mOsm/kg |
| Bacterial screen | PASSED | PASSED |
| Fungal screen | PASSED | PASSED |
| Endotoxin | < 0,05 | < 0,25 EU/mL |
| Mouse Embryo Assay (1-cell) | 100 | ≥ 80 % |

Release date 24 April 2018

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann