

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-1809
Production date	03 December 2018
Expiry date	30 June 2019
Total volume (mL)	41100

	Results	Specifications
pH	7,26	7,20-7,40
Osmolality	281	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date 19 December 2018

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann