

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-1901
Production date	06 February 2019
Expiry date	31 August 2019
Total volume (mL)	42850

	Results	Specifications
pH	7,28	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	85,71	≥ 80 %

Release date 21 February 2019

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann