

Certificate of Analysis GM501 Flush

Batchnumber	FLUSH-2003
Production date	03 April 2020
Expiry date	31 October 2020
Total volume (mL)	27500

	Results	Specifications
pH	7,32	7,20-7,40
Osmolality	279	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date 20 April 2020

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann