

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-2007
Production date	16 September 2020
Expiry date	31 March 2021
Total volume (mL)	67200

	Results	Specifications
pH	7,28	7,20-7,40
Osmolality	283	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	85,71	≥ 80 %

Release date	01 October 2020
--------------	-----------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

p.p. Julia Heinzmann