

Certificate of Analysis GM501 Flush

Batchnumber	FLUSH-2008
Production date	22 October 2020
Expiry date	30 April 2021
Total volume (mL)	58200

	Results	Specifications
pH	7,32	7,20-7,40
Osmolality	283	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Release date	06 November 2020
--------------	------------------

Released by	Dr. Julia Heinzmann
-------------	---------------------

Signature

p.p. Julia Heinzmann