

**Certificate of Analysis
GM501 Flush**

Batchnumber	FLUSH-2102
Production date	17 February 2021
Expiry date	31 August 2021
Total volume (mL)	65600

	Results	Specifications
pH	7,31	7,20-7,40
Osmolality	281	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	81	≥ 80 %

Release date 04 March 2021

Released by Vanessa Bierlich

Signature