

Certificate of Analysis
GM501 Flush

Batchnumber	FLUSH-2202
Production date	09 February 2022
Expiry date	31 August 2022
Total volume (mL)	37500

	Results	Specifications
pH	7,35	7,20-7,40
Osmolality	280	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Release date	24 February 2022
--------------	------------------

Released by	Vanessa Bierlich
-------------	------------------

Signature

P.P. 